

Contents

Part I User Guide

1	Brain–Computer Interfaces	3
1.1	Introduction	3
1.2	Brain–Computer Interfaces (BCIs)	3
1.3	Overview	6
	References	6
2	Brain Sensors and Signals	9
2.1	Relevant Sensors	9
2.2	Brain Signals and Features	10
2.2.1	Using Brain Signals for Communication	10
2.2.2	Mu/Beta Oscillations and Gamma Activity	12
2.2.3	The P300 Evoked Potential	14
2.3	Recording EEG	15
2.3.1	Introduction	15
2.3.2	Electrode Naming and Positioning	15
2.3.3	Important Brain Areas and Landmarks for BCIs	16
2.3.4	Placing Electrodes with a Cap	16
2.3.5	Removing Artifacts and Noise Sources	18
2.4	BCI Signal Processing	21
2.4.1	Introduction	21
2.4.2	Spatial Filtering	22
2.4.3	Feature Extraction: Sensorimotor Rhythms	26
2.4.4	Translation Algorithm	28
	References	29
3	Introducing BCI2000	37
3.1	Motivation	37
3.2	The Design of the BCI2000 Platform	38
3.2.1	A Common Model	38

3.2.2	Source Module and File Format	40
3.2.3	Signal Processing Module	41
3.2.4	User Application Module	41
3.2.5	Operator Module	41
3.2.6	System Variables	42
3.3	BCI2000 Advantages	43
3.4	System Requirements and Real-Time Processing	43
3.5	BCI2000 Implementations and Their Impact	44
	References	45
4	Tour of BCI2000	47
4.1	Starting BCI2000	47
4.1.1	Batch Files	47
4.1.2	BCI2000 Launcher	48
4.2	Getting Help	48
4.3	Configuring BCI2000	49
4.3.1	Module Settings	49
4.3.2	Parameter Files	50
4.3.3	Parameter Help	50
4.4	Important Parameters	51
4.5	Applying Parameters	51
4.6	Starting a Run	52
4.6.1	Data Storage	53
4.7	Retrieving Data Offline	53
4.7.1	Reviewing Brain Signals and Event Markers	53
4.7.2	Viewing Parameters with BCI2000FileInfo	54
4.7.3	Interactions with External Software	55
4.8	Overview of Additional Application Modules	55
4.8.1	Stimulus Presentation	55
4.8.2	P300 Speller	56
4.9	Where to Go from Here	57
5	User Tutorials	59
5.1	General System Configuration	59
5.2	Virtual Cursor Movement with the Sensorimotor Rhythms	59
5.2.1	Obtaining Sensorimotor Rhythm Parameters in an Initial Session	61
5.2.2	Performing the Initial Sensorimotor Rhythm Session	62
5.2.3	Analyzing the Initial Sensorimotor Rhythm Session	63
5.2.4	Configuring Online Feedback	66
5.2.5	Performing a Sensorimotor Rhythm Feedback Session	68
5.3	P300 BCI Tutorial	71
5.3.1	General System Configuration	72
5.3.2	Obtaining P300 Parameters in the Calibration Session	72
5.3.3	Analyzing the Calibration Session with “Offline Analysis”	74

5.3.4	Performing a P300 Spelling Session	79
	References	80
6	Advanced Usage	83
6.1	Matlab MEX Interface	83
6.1.1	Introduction	83
6.1.2	Using BCI2000 MEX Files	83
6.1.3	Building BCI2000 MEX Files	84
6.1.4	BCI2000 MEX Functions	85
6.2	Operator Scripting	87
6.2.1	Events	87
6.2.2	Scripting Commands	88
6.2.3	Examples	88
6.3	Command Line Options	89
6.3.1	Operator Options	89
6.3.2	Core Module Options	89
6.3.3	Data File Formats	90
6.4	AppConnector	92
6.4.1	Introduction	92
6.4.2	Scope	92
6.4.3	Design	93
6.4.4	Description	93
6.4.5	Protocol	94
6.4.6	Examples	94
6.4.7	Parameterization from Within BCI2000	94
6.4.8	Examples	95
6.5	Expression Filter	96
	References	97
7	Programming Reference	99
7.1	Building BCI2000	99
7.1.1	Why Build from Source	99
7.1.2	Tools Required	99
7.1.3	How to Build	100
7.1.4	Starting up BCI2000	101
7.2	Writing a Custom Source Module	101
7.2.1	Example Scenario	102
7.2.2	Writing the ADC Header File	102
7.2.3	ADC Implementation	103
7.2.4	ADC Initialization	104
7.2.5	Data Acquisition	105
7.2.6	Adding the SourceFilter	106
7.2.7	Finished	106
7.3	Writing a Custom Signal Processing Module	106
7.3.1	A Simple Low-Pass Filter	106
7.3.2	The Filter Skeleton	107

7.3.3	The Process Function	107
7.3.4	The Initialize Member Function	108
7.3.5	The Preflight Function	109
7.3.6	Constructor and Destructor	110
7.3.7	Filter Instantiation	111
7.3.8	Visualizing Filter Output	111
7.4	Writing a Custom Matlab Filter	114
7.4.1	Online Algorithm Verification	114
7.4.2	An Example Algorithm in Matlab	114
7.4.3	Transforming Matlab Code into BCI2000 Events	115
8	Exercises	119
8.1	Source Module	119
8.2	Signal Processing Module	121
	References	125
9	Frequently Asked Questions	127
9.1	Timing Issues	127
9.2	Startup	127
9.3	Feedback	128
9.4	Replaying Recorded Data	128
9.5	Random Sequences	129
9.6	Visual Stimulation	129
	References	131

Part II Technical Reference

10	Core Modules	135
10.1	Operator	135
10.1.1	Starting BCI2000	135
10.1.2	Main Window	135
10.1.3	Parameter Configuration Window	138
10.1.4	Matrix Editor Window	140
10.1.5	Visualization Windows	140
10.2	Filter Chain	142
10.3	Data Acquisition Filters	143
10.3.1	DataIOFilter	144
10.3.2	AlignmentFilter	146
10.3.3	TransmissionFilter	147
10.3.4	SourceFilter	147
10.3.5	BCI2000FileWriter	148
10.3.6	EDFFileWriter	149
10.3.7	GDFFileWriter	149
10.3.8	NullFileWriter	151
10.4	Logging Input Filters	151
10.4.1	Input Logger	151

10.4.2	JoystickFilter (Obsolete)	152
10.4.3	KeyLogFilter (Obsolete)	153
10.4.4	MouseFilter (Obsolete)	154
10.5	Signal Source Modules	154
10.5.1	SignalGeneratorADC	154
10.5.2	gUSBampADC	156
10.5.3	gMOBILabADC	159
10.5.4	gMOBILabPlusADC	160
10.6	Signal Processing Filters	162
10.6.1	SpatialFilter	162
10.6.2	ARFilter	165
10.6.3	FFTFilter	166
10.6.4	P3TemporalFilter	168
10.6.5	LinearClassifier	169
10.6.6	Normalizer	172
10.7	Additional Signal Processing Filters	176
10.7.1	LPFilter	176
10.7.2	ConditionalIntegrator	176
10.7.3	StateTransform	177
10.7.4	ExpressionFilter	178
10.7.5	MatlabFilter	179
10.8	Application Modules	180
10.8.1	Cursor Task	180
10.8.2	Stimulus Presentation	184
10.8.3	P300 Speller	191
10.8.4	Keystroke Filter	200
10.8.5	Connector Filters	201
10.9	Tools	201
10.9.1	BCI2000 Offline Analysis	201
10.9.2	USBampGetInfo	205
10.9.3	BCI2000FileInfo	206
10.9.4	BCI2000Export	206
10.10	Localization	208
10.11	P300 Classifier	208
10.11.1	Introduction	208
10.11.2	Interface	209
10.11.3	Reference	213
10.11.4	Tutorial	215
10.11.5	Example	216
	References	217
11	Contributed Modules	219
11.1	Source Modules	219
11.1.1	Amp Server Pro	219
11.1.2	BioRadio	221

11.1.3	BioSemi 2	221
11.1.4	Measurement Computing	222
11.1.5	Data Translation Boards	223
11.1.6	Micromed	225
11.1.7	Modular EEG	226
11.1.8	National Instruments	228
11.1.9	National Instruments MX	228
11.1.10	Neuroscan	229
11.1.11	BrainAmp Systems	231
11.1.12	Tucker-Davis Technologies	232
11.1.13	TMS Refa and Porti Systems	233
11.1.14	BrainProducts V-Amp	234
11.2	Tools	237
11.2.1	EEGLabImport	237
Appendix A The <code>USBampGetInfo</code> Command Line Tool		249
Appendix B The <code>neurogetparams</code> Command Line Tool		253
Index		255

A Practical Guide to Brain-Computer Interfacing with
BCI2000

General-Purpose Software for Brain-Computer Interface
Research, Data Acquisition, Stimulus Presentation, and
Brain Monitoring

Schalk, G.; Mellinger, J.

2010, XXIII, 260 p., Hardcover

ISBN: 978-1-84996-091-5