

EXTRAVASATION KIT

List 4: Contents of extravasation kit

Product	Number	Conditions of storage
Disposable syringes 1 ml 2 ml 5 ml	3 3 3	Room temperature
Disposable cannulas 18 G (pink) 26 G (brown)	5 5	Room temperature
Cold-hot packs in-house product.....	2 for warm compresses 2 for cold compresses	Room temperature Cold: 2–8°C
Gauze pads, sterile, plum size in-house product.....	2 sets at 4 each	Room temperature
Gauze dressings, sterile, approx. 10×10 cm in-house product.....	1 set at 5 pieces	Room temperature
Adhesive plaster in-house product.....	1 roll	Room temperature
Gloves, sterile (medium and large) Protective gloves for handling cytotoxic agents, sterile (medium and large) in-house product.....	1 pair each 1 pair each	Room temperature
Hyaluronidase in-house product.....	1500 IU	Cold: 2–8°C
Dimethylsulfoxide (DMSO) 99% in-house product.....	2 3 100 ml	Room temperature
List 1: Substances List 2: General measures List 3: Substance specific measures List 4: Contents of extravasation kit 1 Extravasation documentation form, 1 tape measure		